

КОДИФИКАТОР**Наборы реагентов**

- | | |
|---|--|
| 251 Вектор-Бест, D-димер-ИФА-Бест | 270 PZ Cormay S.A., D-Dimer kit |
| 271 Abbott, Quantia D-Dimer | 275 Radiometer |
| 253 Alere Triage D-Dimer Test | 283 RAMP D-dimer |
| 273 Beckman Coulter | 261 Roche, Cardiac D-Dimer |
| 255 BioMérieux, VIDAS D-Dimer | 262 Roche, TinaQuant D-Dimer |
| 274 DiaSys, D-Dimer FS | 276 Siemens, Immulite D-Dimer |
| 257 Helena, Auto Blue D-Dimer 400 | 265 Siemens, Innovance D-Dimer |
| 258 Helena, Auto Red D-Dimer 700 | 277 Siemens, Turbo D-Dimer |
| 285 Human, HemoStat D-Dimer | 266 Stago, STA Liatest D-Dimer |
| 259 Instrumentation Laboratory, HemosIL D-Dimer | 268 Tcoag, TriniLIA D-Dimer |
| 260 Instrumentation Laboratory, HemosIL D-Dimer HS | 284 Technoclone, TechnoLEIA D-dimer |
| 286 Mitsubishi, PATHFAST D-Dimer | 279 VEDALAB |
- 99** Другой (опишите в специальном окне электронной формы)

Анализаторы

- | | |
|---|--|
| 1737 Астра, КоаТест-1 | 1708 Olympus AU 400/640/680/2700 |
| 1700 Abbott Architect c4000/c8000/c16000 | 1743 Radim, Alisei Q.S. |
| 1734 Alere, Triage Meter Pro | 1709 Radiometer, AQT 90 Flex |
| 1702 Beckman Coulter AU 400/480/640/2700/5400/5800 | 1723 RAMP Clinical Reader |
| 1719 Biomerieux, Vidas/mini Vidas | 1710 Roche Cardiac Reader |
| 1733 BioSystems, A-15/A-25/BA-400 | 1724 Roche, Cobas C311/C501 |
| 1720 Biotek ELx800/808 | 1729 Roche, Cobas h232 |
| 1703 Clima MC-15 | 1712 Sapphire-400 |
| 1730 Dirui, CS-T240 | 1741 Siemens, BCS XP |
| 1705 Furuno CA-180/270/400 | 1713 Siemens, Immulite 1000/1500/2000 |
| 1706 Helena C-1/C-2/C-4 | 1714 Siemens, SYSMEX CA 500/560/600/1500/7000 |
| 1727 Human, HumaClot Duo/Duo Plus | 1725 Siemens, SYSMEX CS-2000i/2100i/5100i |
| 1707 ILab, ACL Top/Elite/Elite Pro/7000/9000 | 1715 STAGO, STA-R Compact/Evolution/Satellite |
| 1742 Labitec, CoaLAB 1000 | 1735 Stat Fax 2100 |
| 1740 Mindray, MR-96A | 1732 Tecan, Sunrise Absorbance Reader |
| 1726 Mitsubishi, PATHFAST Immunoanalyzer | 1716 Trinity Biotech, AMAX Destiny Plus |
| 1739 Nycocard Reader II | 1717 VEDALAB, Easy Reader |
- 99** Другой (опишите в специальном окне электронной формы)